



T10 FTIR Precision and Bias Review Report

Joe Franklin



- Collected Raw data files from the RR participants
- p Converted the raw files to one software format.
- Performed Spectral Subtraction of New Oil files from Used oil files from Each Lab.
- Analyzed for Oxidation Results based on One macro routine for all files to determine if the variances were purely in the sample prep and data collection.
- p Found that some of the data was more effected by data analysis than others.



- Variance in the data is still coming from the data analysis technique as well as the measurement and sample prep.
- A full workshop is required to improve the precision of the test and it can be done.
- p I will set up this workshop in the next 30 days.