



## Test Monitoring Center

6555 Penn Avenue  
Pittsburgh, PA 15206-4489  
(412) 365-1000

MEMORANDUM: 00-150  
DATE: October 23, 2000  
TO: D5800 (EVLO) Participants  
FROM: Tom Schofield  
SUBJECT: D5800 Test Method Update

I have been advised that a revised D5800 test method is available for purchase from ASTM Headquarters. The newest update is D5800-00a, and contains all three Noack procedures (Woods Metal, Non-Woods Metal & Selby Noack). Effective immediately, the TMC will be monitoring the D5800-00a for future calibrations.

Please be aware that, of this writing, the ASTM website ([www.astm.org](http://www.astm.org)) still lists D5800-00 as the latest available version. D5800-00 is NOT the same as D5800-00a. The version with the lower-case "a" contains all three procedures while D5800-00 does not. Please be explicit when ordering the new method. Until the ASTM website is updated, D5800-00a is available by contacting ASTM Headquarters by phone at 610-832-9500. You could also contact them by fax (610-832-9555) or e-mail ([service@astm.org](mailto:service@astm.org)).

Please contact me if you require additional information.

TMS/tms

Distributed via Internet mail message

c: <ftp://tmc.astm.cmri.cmu.edu/docs/bench/d5800/memos/mem00-150>